



Determination of Four Chemical Characterisation Compounds in Honey by Liquid Chromatography Tandem Mass Spectrometry (LC- MS/MS)

Chemistry Laboratory Method

MPI Technical – Paper No: 2017/30

ISBN No: 978-1-77665-435-2 (online)

ISSN No: 2253-3923 (online)

April 2017

CLM-HON1.09	Page 2 of 30
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS	
Revision: 09	Replaces: CLM-HON1.08 Effective: 10/04/2017

Disclaimer

While every effort has been made to ensure the information in this publication is accurate, the Ministry for Primary Industries does not accept any responsibility or liability for error of fact, omission, interpretation or opinion that may be present, nor for the consequences of any decisions based on this information.

Requests for further copies should be directed to:

Publications Logistics Officer
Ministry for Primary Industries
PO Box 2526
WELLINGTON 6140

Email: brand@mpi.govt.nz
Telephone: 0800 00 83 33
Facsimile: 04-894 0300

This publication is also available on the Ministry for Primary Industries website at <http://www.mpi.govt.nz/news-and-resources/publications/>

© Crown Copyright - Ministry for Primary Industries

CLM-HON1.09	Page 3 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

Contents		Page
1	Introduction	5
1.1	Background	5
1.2	Summary of Procedure	5
1.3	Applicability	5
2	Equipment	5
2.1	Apparatus	5
2.2	Instrumentation	6
3	Reagents and Solutions	6
3.1	Reagents	6
3.2	Solutions	6
4	Standard(s)	7
4.1	Standard Information	7
4.2	Preparation of Standard Solution(s)	7
5	Sample Preparation	8
5.1	Sample Homogenisation	8
6	Analytical Procedure	8
6.1	Preparation of Controls	8
6.2	Sample Extraction	9
6.3	Preparation of Calibration Curve Standards	9
6.4	Instrument Settings	10
6.5	Injection Sequence/Sample Set	12
7	Calculations / Identification	13
7.1	Calculations	13
7.2	Quantification Calculation	13
7.3	Confirmation Criteria	14
7.4	Quantification Criteria	14
8	Safety Information and Precautions	14
8.1	Required Protective Equipment	14
8.2	Hazards	14
9	Quality Assurance Plan	14
9.1	Performance Standard	14
9.2	Critical Control Points and Specifications	15
9.3	Inter-laboratory Check Samples	15
9.4	Condition Upon Receipt	15
9.5	Reporting	15
10	Appendix	16
10.1	Sample Chromatograms and Calibration Curves	16

CLM-HON1.09	Page 4 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

10.2 Structures	27
10.3 Validation Data	28
10.4 Synthetic Honey Substitute for Controls	29
11 Approvals and Authorities	30
11.1 Approvals on File	30

CLM-HON1.09	Page 5 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

1 Introduction

1.1 BACKGROUND

This method was developed and validated by the National Measurement Institute, Port Melbourne, Australia as part of the Ministry for Primary Industries Mānuka Honey Science programme. The assay is one of the key tests that will be used to verify authenticity of monofloral and multifloral Mānuka honey.

The validated assay was used to determine the specificity and concentrations of the target chemicals in nectar (~500 samples) from a variety of plant species representing 17 *Leptospermum* species, 5 *Kunzea* species and 13 other plant species relevant to honey production in New Zealand. The assay has been tested on over 800 honey samples from various floral sources and geographic locations. Samples were sourced from two New Zealand flowering seasons, New Zealand industry archives (up to 5 years) and foreign sourced honeys (e.g. Australia, China, Africa, USA, Europe).

1.2 SUMMARY OF PROCEDURE

The four chemical characterisation compounds are extracted from honey with a mixture of aqueous acetonitrile and formic acid. The extract is further diluted and filtered before determination by LC-MS/MS.

1.3 APPLICABILITY

This method is suitable for the quantification of chemical characterisation compounds in honeys of various floral types at the levels listed in Table 1.

Table 1: Chemical Characterisation Compounds

Compound Common Name	Abbreviation	CAS#	Limit of Reporting (LOR) (mg/kg)	Linear Range (mg/kg)
2'-Methoxyacetophenone	2MAP	579-74-8	1	1 – 100
2-Methoxybenzoic acid	2MBA	579-75-9	1	1 – 100
3-Phenyllactic acid	3PA	828-01-3	1	1 – 100
4-Hydroxyphenyllactic acid	4HPA	306-23-0 or 6482-98-0	1	1 – 100
Forchlorfenuron (<i>internal standard only</i>)	FCF	68157-60-8	-	-

2 Equipment

Note: Equivalent equipment may be substituted.

2.1 APPARATUS

- Balance – analytical ± 0.0001 g
- Balance – top loading ± 0.01 g
- Bottle – amber, glass, 10 mL

CLM-HON1.09	Page 6 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

- d. Beakers – 10 mL
- e. Centrifuge – refrigerated, with rotor for 50 mL tubes
- f. Centrifuge tubes – 15 and 50 mL conical, disposable, polypropylene with caps
- g. Pipettes – Adjustable 10 – 5000 µL
- h. Shaker – horizontal flatbed
- i. Syringe filter, 13 mm diameter 0.45 µm PTFE
- j. Vials – amber, glass 1.5 or well plate suitable for use in an auto-sampler
- k. Volumetric flasks – 10, 50, 100 mL and 1 L
- l. Vortex mixer – variable speed

2.2 INSTRUMENTATION

- a. Waters Acquity UPLC I-Class System
- b. Waters Xevo TQ-S Triple Quadrupole MS
- c. HPLC Guard Column – Waters Acquity UPLC BEH C18 VanGuard Pre-column, 130Å, 1.7 mm, 2.1 mm x 5 mm, (Waters P/N 186003975)
- d. HPLC Column – Waters Acquity UPLC BEH C18 Column, 130Å, 1.7 µm, 2.1 mm X 100 mm, (Waters P/N 186002352)

3 Reagents and Solutions

Note: All reagents are AR Grade, unless otherwise specified. Equivalent reagents / solutions may be substituted. The stability time frame of the solution is dependent on the expiration dates of the compounds used. The maximum length of time that a working reagent shall be used is 1 year unless the laboratory has produced extension data.

3.1 REAGENTS

- a. Acetonitrile – HPLC grade
- b. Formic acid – Fisher Chemical, 99.5+%, Optima™ LC/MS Grade, P/N A117-50
- c. Water – Deionized, HPLC grade

3.2 SOLUTIONS

a. **Extraction solution: 10:1:90 Acetonitrile/formic acid/water**

Approximately half fill a 1L volumetric flask with deionised water. Measure 100 mL of acetonitrile and measure 10 mL of formic acid into the volumetric flask. Fill to volume with deionised water and invert gently to mix.

Store at room temperature. Prepare fresh weekly.

b. **1% Formic acid in acetonitrile**

Approximately half fill a 100 mL volumetric flask with acetonitrile. Pipette 1 mL of formic acid into the volumetric flask. Fill to volume with acetonitrile and invert gently to mix.

Store at room temperature. Prepare fresh for each use.

CLM-HON1.09	Page 7 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

c. Mobile Phase A: 0.1% Formic acid in water

Approximately half fill a 1L volumetric flask with deionised water. Pipette 1 mL of formic acid into the volumetric flask. Fill to volume with deionised water and invert gently to mix.

Store at room temperature. Prepare fresh for each analytical batch.

d. Mobile Phase B: 0.1% Formic acid in acetonitrile

Approximately half fill a 1L volumetric flask with acetonitrile. Pipette 1 mL of formic acid into the volumetric flask. Fill to volume with acetonitrile and invert gently to mix.

Store at room temperature. Prepare fresh for each analytical batch.

4 Standard(s)

Note: Equivalent standards / solutions may be substituted. Purity and counter-ions are to be taken into account when calculating standard concentrations. The stability time frame of the solution is dependent on the expiration date of the components used. In-house prepared standards shall be assigned an expiration date that is no later than the expiration date of the earliest expiring component or no later than the stability stated in the method, whichever ends soonest. The maximum length of time that an in-house prepared standard shall be used is 3 months unless the laboratory has produced extension data.

4.1 STANDARD INFORMATION

- 2'-Methoxyacetophenone – 99% pure. Catalogue No. M9203 Sigma-Aldrich
- 2-Methoxybenzoic acid – 99% pure. Catalogue No. 169978 Sigma-Aldrich
- 3-Phenyllactic acid – 97% pure. Catalogue No. P7251 Sigma-Aldrich
- 4-Hydroxyphenyllactic acid – 95% pure. Catalogue No. 222023 Fluorochem
- Forchlorfenuron (Internal Standard) – 99% pure. Catalogue No. 32974 Sigma-Aldrich

4.2 PREPARATION OF STANDARD SOLUTION(S)

Note: Adjust all standard weights for purity.

a. Stock standard and internal standard solutions (~1000/10000 mg/L)

Weigh approximately 10 mg of the 2MAP, 2MBA, 4HPA and FCF standards into its own 10 mL volumetric flask and bring to volume with acetonitrile. Weigh approximately 100 mg of the 3PA standard into its own 10 mL volumetric flask and bring to volume with acetonitrile. Record the weight to 0.1 mg and calculate the exact concentration, taking into account both the moisture, purity and salt content. Transfer to an amber, glass bottle with minimal headspace.

These standards are stable for 3 months when stored in a refrigerator at 2 – 8°C.

CLM-HON1.09	Page 8 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

b. Mixed spiking standard solution (100/1000 mg/L)

Pipette ~1000 μL (adjusted for the actual stock standard concentration) of each stock standard solution into a 10 mL volumetric flask and bring to volume with 1% formic acid in acetonitrile. Transfer this solution into a 10 mL amber, glass bottle with minimal headspace.

This standard is stable for 3 months when stored in a refrigerator at 2 – 8°C.

c. Mixed spiking standard solution (10/100 mg/L)

Pipette 1000 μL of the **100/1000 mg/L** mixed spiking standard solution into a 10 mL volumetric flask and bring to volume with 1% formic acid in acetonitrile. Transfer this solution into a 10 mL amber, glass bottle with minimal headspace.

This standard is stable for 3 months when stored in a refrigerator at 2 – 8°C.

d. Mixed calibration standard solution (0.25 and 2.5 mg/L)

Pipette 250 μL of the **10/100 mg/L** mixed spiking standard solution into a 10 mL volumetric flask and bring to volume with 1% formic acid in acetonitrile. Transfer this solution into a 10 mL amber, glass bottle with minimal headspace.

This standard is stable for 3 months when stored in a refrigerator at 2 – 8°C.

e. Forchlorfenuron internal standard solution (~10 mg/L)

Pipette 500 μL of the stock FCF internal standard solution into a 50 mL volumetric flask and bring to volume with 1% formic acid in acetonitrile. Transfer this solution into a 10 mL amber, glass bottle with minimal headspace.

This standard is stable for 3 months when stored in a refrigerator at 2 – 8°C.

5 Sample Preparation

5.1 SAMPLE HOMOGENISATION

- Ensure honey samples are brought to room temperature and thoroughly mixed before sub-sampling.

6 Analytical Procedure

6.1 PREPARATION OF CONTROLS

- Weigh four $1\text{ g} \pm 0.1\text{ g}$ of blank homogenised honey portions into 50 mL disposable centrifuge tubes.

Note: A honey known to contain negligible concentrations of the four characterisation compounds should be used for the control samples.

- Fortify the positive control samples with the volumes of mixed spiking solution as described in Table 2.

CLM-HON1.09	Page 9 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

Table 2: Fortification of positive control samples

Positive Control and Concentration	Concentration of Mixed Spiking Standard Solution	Volume of Mixed Spiking Standard Solution (μL)
LOR ¹ (1/10 mg/kg)	10/100	100
2x LOR (2/20 mg/kg)	10/100	200
10x LOR (10/100 mg/kg)	100/1000	100

6.2 SAMPLE EXTRACTION

- a. Weigh 1 g \pm 0.1 g of homogenised honey samples into 50 mL disposable centrifuge tubes. Record the weight of the analytical portion to 0.01 g accuracy. Include the prepared controls in the sample set at this time.

Note: In addition to the control samples, one of the analytical samples is weighed in duplicate per batch. Also it is recommended laboratories develop and run an in-house quality control (QC) sample for batch to batch perform monitoring.

- b. Add 9 mL of 10:1:90 acetonitrile/formic acid/water solution to each tube.
c. Vortex mix vigorously to ensure thorough mixing.
d. Agitate the tubes on a horizontal shaker vigorously for 10 minutes.
e. Centrifuge the tubes for 10 minutes at approximately 3500 RCF at 10°C.
f. Filter ~ 1 mL of each dilution through a 0.45 μm PTFE filter into a 1.5 mL, amber auto-sampler vial.

6.2.1 200 times dilution (controls and samples)

- g. Add 500 μL of the supernatant (step 6.2 f.), 100 μL of forchlorfenuron internal standard solution (10 mg/L) and 9.4 mL of 10:1:90 acetonitrile/formic acid/water solution into a 15 mL, polypropylene tube. Cap and mix thoroughly.

6.2.2 5000 times dilution (samples only)

- h. Add 20 μL of the supernatant (step 6.2 f.), 100 μL of forchlorfenuron internal standard solution (10 mg/L) and 9.88 mL of 10:1:90 acetonitrile/formic acid/water solution into a 15 mL, polypropylene tube. Cap and mix thoroughly.
i. Transfer ~ 1 mL of each dilution into a 1.5 mL, amber auto-sampler vial for LC-MS/MS determination

Note 1: The dilution series above are recommended based on the instrument sensitivity achieved at the time of validation. Laboratory may need to adjust dilution ratios depending on instrument performance.

Note 2: Typically, both recommended dilution series are required to be analysed to cover the expected range of characterisation compounds found in Mānuka honey.

Note 3: In addition, a reagent blank of acetonitrile solution is prepared and injected the instrument sequence between the positive controls and the sample set.

6.3 PREPARATION OF CALIBRATION CURVE STANDARDS

- a. Prepare the calibration standards using 10:1:90 acetonitrile/formic acid/water solution according to table 3 in 1.5 mL amber, glass vials. Cap and vortex mix. Prepare calibration curve standards fresh for every batch and store at 2 – 8°C.

¹ Note: The LOR referred to corresponds to 2MAP, 2MBA, and 4HPA. The spiking level for 3PA is 10 times that compound's LOR.

CLM-HON1.09	Page 10 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

Table 3: Preparation of calibration curve standards

2MAP, 2MBA, 4HPA Cal. Std. Conc. (mg/L)	3PA Cal. Std. Conc. (mg/L)	2MAP, 2MBA, 4HPA Honey Conc. (mg/kg)	3PA Honey Conc. (mg/kg)	Volume of Mixed Cal. Standard (μ L)	Volume of FCF Internal Standard (μ L)	Volume. of Extraction Solution (μ L)
0	0	0	0	0	10	990
0.0025	0.025	0.5	5	10	10	980
0.005	0.05	1	10	20	10	970
0.01	0.1	2	20	40	10	950
0.025	0.25	5	50	100	10	890
0.05	0.5	10	100	200	10	790
0.1	1	20	200	400	10	590

6.4 INSTRUMENT SETTINGS

Note: The instrument parameters may be optimized to ensure system suitability.

Table 4: HPLC Conditions:

Mobile Phase A	0.1% Formic acid in water
Mobile Phase B	0.1% Formic acid in acetonitrile
Flow Rate	0.2 mL/min
Column Temperature	40°C \pm 2°C
Injection Volume	2 μ L
Auto-sampler Temperature	8°C \pm 5°C
Run Time	10 minutes

Table 5: HPLC Run Events – Diverter value

Event Sequence	Event Description	Time (mins)
Event 1	Flow State Waste	0.00
Event 2	Flow State LC	3.20
Event 3	Flow State Waste	5.95

Table 6: HPLC Mobile Phase Gradient Table:

Time (mins)	% Mobile Phase A	% Mobile Phase B
0.00	95	5
0.75	95	5
2.00	85	15
4.00	30	70
6.00	2	98
6.50	2	98
7.00	95	5
10.00	95	5

CLM-HON1.09	Page 11 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

Table 7: Instrument Set-Up and Parameters

Negative Ion Mode (ES-) 4-Hydroxyphenyllactic acid 3-Phenyllactic acid Forchlorfenuron (internal standard)		Positive Ion Mode (ES+) 2'-Methoxyacetophenone Methoxybenzoic acid	
Capillary voltage	2.00 kV	Capillary voltage	3.00 kV
Cone voltage	40.00 V	Cone voltage	40.00 V
Source Temperature	150°C	Source Temperature	150°C
Desolvation Temperature	450°C	Desolvation Temperature	450°C
Cone Gas Flow	150 L/Hr	Cone Gas Flow	150 L/Hr
Desolvation Gas Flow	800 L/Hr	Desolvation Gas Flow	800 L/Hr
Collision Gas Flow	0.14 mL/min	Collision Gas Flow	0.13 mL/min
Nebuliser Gas Flow (Bar)	7	Nebuliser Gas Flow (Bar)	7
LM 1 Resolution	2.9	LM 1 Resolution	2.8
HM 1 Resolution	14.9	HM 1 Resolution	14.9
Ion Energy 1	0.4	Ion Energy 1	0.5
MS Mode Entrance	1	MS Mode Entrance	1
MS Mode Collision Energy	5	MS Mode Collision Energy	5
MSMS Mode Collision Energy	20	MSMS Mode Collision Energy	20
MS Mode Exit	1	MS Mode Exit	1
MSMS Mode Entrance	30	MSMS Mode Entrance	30
MSMS Mode Collision Energy	20	MSMS Mode Collision Energy	20
MSMS Mode Exit	30	MSMS Mode Exit	30
LM 2 Resolution	2.8	LM 2 Resolution	2.7
HM 2 Resolution	15	HM 2 Resolution	15
Ion Energy 2	0.7	Ion Energy 2	0.5
Gain	1	Gain	1

Table 8: MS/MS Parameters

Compound	Retention window (min.)	Precursor Ion [†] (m/z)	Product Ion [†] (m/z)*	Accurate Mass Observed ‡	Dwell Time [†] (Auto) (sec)	Cone Voltage [†] (Volts)	Collision Energy [†] (Volts)
2MAP	5.3 – 6.0	150.70	79.31	151.0752	0.029	2	24
		150.70	105.32	79.0546 105.0700			
2MBA	4.7 – 5.3	152.99	77.07	153.0545	0.019	30	39
		152.99	92.02	77.0389 92.0259			
3PA	4.6 – 5.0	165.01	103.13	165.0548	0.019	45	18
		165.01	119.12	103.0543 119.0492			
4HPA	3.8 – 4.2	181.01	73.11	181.0498	0.025	50	14
		181.01	119.12	72.9920 119.0492			
		181.01	135.05	135.0441			
FCF (Internal Standard)	5.4 – 5.9	245.82	90.95	246.0441	0.036	4	28
		245.82	126.87	91.0294 127.0063			

* Most abundant product ion (quantification ion) is in bold.

† The parameters shown in the table are specific to instrumentation used at the time of validation. These can vary depending on the instrumentation used and the most appropriate precursor and product ions needs to be determined by each laboratory.

‡ The observed accurate masses shown in the table were found using liquid chromatography high-resolution accurate-mass spectrometry.

6.5 INJECTION SEQUENCE/SAMPLE SET

- Calibration Curve
- Positive Controls
- Matrix Blank
- Reagent Blank
- QC sample
- Samples and duplicates

Calibration standards should be injected at the start of each batch and one of these standards re-injected after every 8 – 10 samples throughout the batch.

The lowest calibration standard point (0.0025/0.025 mg/L) is to be re-injected at the end of each analytical batch sequence as a book end to enable the analyst to check the consistency and stability of the instrument's response. If instrument response is found to have significantly degraded (or enhanced) then the analyst must consider the impact to the veracity of the findings and investigate potential causes.

Linearity of standards has been checked and found to be generally good in the range equivalent to sample concentrations 1 – 100 mg/kg. This can vary depending on instrumentation used and needs to be confirmed by each laboratory. Samples with

CLM-HON1.09	Page 13 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

calculated concentration outside the linear range of the methods should be diluted to a suitable concentration.

7 Calculations / Identification

7.1 CALCULATIONS

a. Calculation of batch recovery

$$\text{Recovery (\%)}_{\text{Spike}} = \frac{\text{Conc. Spike} - \text{Conc. Blank}}{\text{Conc. Added}} \times 100$$

where:

Recovery (%) Spike = percentage recovery of the control spike

Conc. Spike = calculated concentration of the control spike

Conc. Blank = calculated concentration of the matrix blank

Conc. Added = expected concentration of the control spike

The percentage recovery is found for each of the each spiking levels

b. Calculation of relative percentage difference (RPD) of duplicates

$$\text{RPD (\%)} = \frac{|\text{Conc. Duplicate 1} - \text{Conc. Duplicate 2}|}{\text{Mean Conc. (Duplicate 1 \& 2)}} \times 100$$

where:

RPD (%) = relative percentage difference of the duplicates

Conc. Duplicate 1 = calculated concentration of the first duplicate

Conc. Duplicate 2 = calculated concentration of the second duplicate

Mean Conc. (Duplicate 1 & 2) = arithmetic mean of the calculated concentration of the duplicates

7.2 QUANTIFICATION CALCULATION

- Peak areas of analytes and internal standards are used for quantification.
- The coefficient of correlation (r^2) must be ≥ 0.99 .
- Do not use the origin as a data regression point.
- Determine sample concentrations for 2MAP, 3PA and 4HPA is typically found using a weighted linear regression ($1/x$) and determine sample concentrations for 2MBA is typically found using a weighted quadratic regression ($1/x$). This can vary depending on the instrumentation used and the most appropriate curve fitting needs to be determined by each laboratory. It is recommended curve fitting is optimised to minimise residuals for each calibration point.
- Calculation of concentration of characterisation compounds is found as follows.
(linear regression example)

$$\text{Calculated concentration (mg/kg)} = \frac{(y - c)}{m} \times d \times \frac{1}{S_w}$$

where:

y = ratio of the quantification ion peak areas of the analyte/internal standard

CLM-HON1.09	Page 14 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

c = y intercept
 m = gradient
 d = dilution factor (typically 200 or 5000)
 s_w = sample weight (g)

7.3 CONFIRMATION CRITERIA

- The retention time of the analyte(s) must match that of the positive control or the external standard injected most recently before the relevant sample within ± 0.1 min.
- Product ion abundance ratios must match that of the positive control or the external standard injected most recently before the relevant sample within $\pm 30\%$ (relative).
- Each ion must have a signal-to-noise ratio ≥ 3 .

7.4 QUANTIFICATION CRITERIA

- The sample peak retention time must be within ± 0.1 min of the positive control or the external standard injected most recently before the relevant sample.
- The quantification ion must have a signal to noise ratio of ≥ 10 .
- The additional ion(s) must be present in sample with a signal to noise ratio of ≥ 3 .
- The coefficient of determination (r^2) for the calibration curve must be ≥ 0.99

Note: Quantification criteria are required only for analytes that are to be quantified in the sample set.

8 Safety Information and Precautions

8.1 REQUIRED PROTECTIVE EQUIPMENT

- Safety glasses and/or face shield
- Disposable gloves
- Laboratory coat

8.2 HAZARDS

Table 9: Hazards

Procedure Step	Hazard	Recommended Safe Procedures
Acetonitrile	Flammable and poisonous	Use reagents in an efficient fume hood away from all electrical devices and open flames. Wear gloves and protective eyewear.
Formic acid	Acid burns	Wear protective equipment and avoid contact with skin.

9 Quality Assurance Plan

9.1 PERFORMANCE STANDARD

- Positive controls are positive for all analytes using the criteria in Section 7.
- The positive control recovery at the LOR¹ spiking level is between 80 – 110% for all analytes that will be quantified. For any recovery $<80\%$ or $>110\%$, the data should be reviewed and/or the batch repeated.

CLM-HON1.09	Page 15 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

- c. The calculated amounts for each of the compounds of the QC sample should be in an acceptable range (mean \pm 2 standard deviations) determined by using a control chart.
- d. The relative percent difference (RPD) between duplicates should be <20%.
- e. Relative response factors for standards injected through the analytical run should vary by no more than 15%.

9.2 CRITICAL CONTROL POINTS AND SPECIFICATIONS

<u>Record</u>	<u>Acceptable Control</u>
Sample weight of honey	1.0 ± 0.1 g

9.3 INTER-LABORATORY CHECK SAMPLES

- a. System, minimum contents.
 - i. Frequency: As available, per analyst when samples tested.
 - ii. Records are to be maintained.
- b. Acceptability criteria.
Refer to 9.1.
If unacceptable values are obtained, then:
 - i. Investigate following established procedures.
 - ii. Take corrective action as warranted.

9.4 CONDITION UPON RECEIPT

- a. Room temperature, no evidence of spoilage, leakage or container damage.

9.5 REPORTING

- a. Each characterisation compound is reported in units of mg/kg to two significant figures.

10 Appendix

10.1 SAMPLE CHROMATOGRAMS AND CALIBRATION CURVES

Following are a series of chromatograms obtained for analysis of honey prepared and analysed during the initial validation of the method.

Figure 1: Blank in positive and negative ion






Figure 2: 0.002 mg/L Matrix standards

0.002ppm Matrix standard


MH14_161024_006 Sm (Mn, 3x2)




MH14_161024_006 Sm (Mn, 3x2)



MH14_161024_006 Sm (Mn, 3x2)



MH14_161024_006 Sm (Mn, 3x2)



MH14_161024_006 Sm (Mn, 3x2)




Figure 3: 1.0 mg/L Matrix standards**1.0ppm Matrix standard**

MH14_161024_011 Sm (Mn, 3x2)




Figure 4: 1 mg/kg Honey spike

1.0ppm Spike

MH14_161024_015 Sm (Mn, 3x2)




Figure 5: 100 mg/kg Honey spike

100ppm Spike

MH14_161024_019 Sm(Mn, 3x2)




Figure 6: MH2 - Commercial source honey in negative ion (TIC and MRM)





Figure 7: MRM of 2 Channels ES- for 3-phenyllactic acid of 1.0 mg/kg honey spike and 0.05 mg/kg matrix standard

Figure 8: MH2- Commercial source honey in positive ion (TIC and 3MRM ES+ 2-methoxybenzoic acid)








Figure 9: MRM of 3 Channels ES+ for 2'-Methoxyacetophenone of 1.0mg/kg honey spike and 1.0 mg/kg matrix standard

Figure10: Standard calibration curves for four chemical characterisation compounds of honey in the range of 0.002 mg/L - 1 mg/L.

Compound name: DL-3-Phenylactic Acid
Correlation coefficient: $r = 0.999919$, $r^2 = 0.999838$
Calibration curve: $4.4103 * x + 0.00208755$
Response type: Internal Std (Ref7), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None







Compound name: 4-HydroxyPhenylactic Acid
Correlation coefficient: $r = 0.999693$, $r^2 = 0.999386$
Calibration curve: $3.10816 * x + 0.00147736$
Response type: Internal Std (Ref7), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



CLM-HON1.09	Page 27 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

10.2 STRUCTURES

Table 10: Characterisation compounds detail

Compound	CAS No.	Molecular Formula	Structure
2'-Methoxyacetophenone	579-74-8	C ₉ H ₁₀ O ₂	
2-Methoxybenzoic acid	579-75-9	C ₈ H ₈ O ₃	
3-Phenyllactic acid	828-01-3	C ₉ H ₁₀ O ₃	
4-Hydroxyphenyllactic acid	306-23-0 / 6482-98-0	C ₉ H ₁₀ O ₄	

CLM-HON1.09	Page 28 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

10.3 VALIDATION DATA

- a. Spiked recoveries at LOR, 2xLOR, 5xLOR, 50xLOR, 100xLOR, were carried out between three technicians from 01/08/2016 – 25/08/2016.

Table 11: Mean recoveries of honey spiked with the four characterisation compounds (1 mg/kg – 100 mg/kg)

Compound	n	Recovery (%)	Standard Deviation (SD) (%)	Coefficient of Variation (CV) (%)	Acceptable Range ($\bar{x} \pm 2SD$) (%)
2'-Methoxyacetophenone	40	102	7.9	7.8	86 – 117
2-Methoxybenzoic acid	40	102	6.8	6.7	88 – 115
3-Phenyllactic acid	40	90	6.9	7.7	76 – 104
4-Hydroxyphenyllactic acid	40	91	7.2	7.9	77 – 105

Table 12: Summary of individual recoveries at all levels for honey spiked with the four characterisation compounds

Compound	Spike Level (mg/kg)	n	Recovery (%)	Range (%)	SD (%)	CV (%)
2'-Methoxyacetophenone	1	8	98	86 – 110	5.9	6.0
	2	8	94	84 – 104	8.1	5.4
	5	8	107	92 – 121	7.3	6.9
	50	8	108	92 – 123	7.9	7.3
	100	8	102	94 – 111	4.2	4.1
2-Methoxybenzoic acid	1	8	106	94 – 118	5.9	5.6
	2	8	105	89 – 120	7.9	7.5
	5	8	102	92 – 112	5.1	5.0
	50	8	99	89 – 108	4.7	4.7
	100	8	98	83 – 112	7.2	7.4
3-Phenyllactic acid	1	8	90	79 – 101	5.3	5.8
	2	8	83	75 – 91	3.8	4.6
	5	8	91	80 – 102	5.3	5.8
	50	8	93	78 – 108	7.5	8.1
	100	8	92	75 – 109	8.6	9.4
4-Hydroxyphenyllactic acid	1	8	97	82 – 112	7.6	7.9
	2	8	87	78 – 96	4.3	5.0
	5	8	92	79 – 105	6.3	6.9
	50	8	90	76 – 104	7.2	8.0
	100	8	90	74 – 106	7.9	8.8

- b. Measurement uncertainty determination method

The combined standard uncertainty was determined by considering the following major sources of measurement uncertainty; sample recovery, sample homogeneity, purity and linearity of analytical standard. A coverage factor of 2 (K=2) was applied to all combined

CLM-HON1.09	Page 29 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

standard uncertainties to give a 95% confidence interval (CI) for the % measurement uncertainty estimates.

Homogeneity: Relative Standard Deviation (RSD) of duplicate analyses of 22 honey samples for each of the characterisation compounds.

Table 14: Relative standard deviation (RSD) for homogeneity estimates

Compound	n	RSD
2'-Methoxyacetophenone	22	0.0177
2-Methoxybenzoic acid	22	0.0184
3-Phenyllactic acid	22	0.0185
4-Hydroxyphenyllactic acid	22	0.0105

Method Recovery: The RSD of all recoveries from 3 batches carried out by three different analysts for concentrations from 1 to 100 mg/kg.

Calibration: The RSD of the calibration curve (NMI calculation template). This may be a slight overestimation as it has been determined using non-weighted concentrations whereas the calibration used by the instrument software has a $\frac{1}{x}$ weighting.

Standard Mass: The RSD of the standard mass for preparation of the stock standards.

Other sources: The total from all other sources estimated using scientific judgement.

Results are corrected for recovery. Therefore a normal distribution can be assumed for the remaining combined uncertainties.

c. The measurement uncertainty was estimated to be as follows for honey matrices. The resulting combined uncertainties are displayed in the tables below:

Table 14: Measurement uncertainties

Compound	LOD (mg/kg)	LOR (mg/kg)	Linear Range (mg/kg)	Measurement uncertainty at the 95% confidence level and concentration at which MU was established
2'-Methoxyacetophenone	0.5	1	1 – 100	± 20% at 5 mg/kg
2-Methoxybenzoic acid	0.5	1	1 – 100	± 25% at 5 mg/kg
3-Phenyllactic acid	0.5	1	1 – 100	± 20% at 5 mg/kg
4-Hydroxyphenyllactic acid	0.5	1	1 – 100	± 20% at 5 mg/kg

10.4 SYNTHETIC HONEY SUBSTITUTE FOR CONTROLS

If no naturally blank honey is available, a synthetic honey substitute can be used as an alternative blank matrix.

Note: The use of nature honey as blank matrix is recommended where possible.

CLM-HON1.09	Page 30 of 30	
Title: Determination of Four Chemical Characterisation Compounds in Honey by LC-MS/MS		
Revision: 09	Replaces: CLM-HON1.08	Effective: 10/04/2017

10.4.1 Equipment

- a. Beaker, glass, 50 mL
- b. Hotplate, stirrer

10.4.2 Reagents

- a. D-(-)-Fructose, $\geq 99\%$ – ACS reagent grade, Sigma-Aldrich F0127
- b. D-(+)-Glucose, $\geq 99.5\%$ – ACS reagent grade, Sigma-Aldrich G8270
- c. D-(+)-Maltose, $\geq 99\%$ – ACS reagent grade, Sigma-Aldrich M5885
- d. Sucrose, $\geq 99.5\%$ – ACS reagent grade, Sigma-Aldrich S0389
- e. Water, deionised

10.4.3 Method

Into a 50 mL beaker weigh $8.4 \text{ g} \pm 0.1 \text{ g}$ of fructose, $6.6 \text{ g} \pm 0.1 \text{ g}$ of glucose, $0.26 \text{ g} \pm 0.01 \text{ g}$ of maltose, $0.12 \text{ g} \pm 0.01 \text{ g}$ of sucrose and $4.6 \text{ g} \pm 0.1 \text{ g}$ of deionised water. Place on a hotplate stirrer and warm gently while mixing until all sugars are completely dissolved. Remove from heat and allow to cool. Transfer to an air-tight container. Stable for 6 months when stored in a refrigerator at $2 - 8^\circ\text{C}$.

11 Approvals and Authorities

11.1 APPROVALS ON FILE

- a. Issuing Authority: Manager, Chemical and Microbiological Assurance